

KRONOS 18 L

RECOMMENDED PROGRAMS

● OPEN ● WRAPPED

B-CLASS

TEMPERATURE AND NAME OF CYCLE	TYPE	TOTAL TIME	WEIGHT	STERILIZABLE MATERIALS AND INSTRUMENTS	WRPD
HELIX TEST - B&D TEST	Test	22 M	0	Steam penetration test	-
VACUUM TEST		15 M		Leak detection test	-
121° C HOLLOW WRAPPED	B	50 M	2,5	Delicate hollow, stainless steel hollow and turbines	●
134° C HOLLOW WRAPPED		39 M		Stainless steel hollow and turbines	●
134° C PRION		54 M		Stainless steel hollow and metal instruments	●
134° C POROUS		50 M	1,0	Porous loads	●
134° C HOLLOW OPEN		36 M	4,0	Stainless steel hollow and turbines	●

REFERRING TO KRONOS B18

S-CLASS

TEMPERATURE AND NAME OF CYCLE	TYPE	TOTAL TIME	WEIGHT	STERILIZABLE MATERIALS AND INSTRUMENTS	WRPD
BOWIE & DICK TEST	Test	28 M	0	Steam penetration test	-
VACUUM TEST		15 M		Leak detection test	-
121° C HOLLOW WRAPPED	S	56 M	1,5	Delicate hollow, stainless steel hollow	●
134° C HOLLOW WRAPPED		45 M		Stainless steel hollow	●
134° C PRION		60 M		Stainless steel hollow and metal instruments	●
134° C POROUS		56 M	0,5	Little porous loads	●
134° C FAST		36 M	2,5	Rubber solids and metal solids	●

REFERRING TO KRONOS S18

N-CLASS

TEMPERATURE AND NAME OF CYCLE	TYPE	TOTAL TIME	WEIGHT	STERILIZABLE MATERIALS AND INSTRUMENTS	WRPD
121° C SOLIDS OPEN	N	48 M	1,5	Rubber solids and delicate solids	●
134° C SOLIDS OPEN		36 M		Rubber solids and metal solids	●
134° C PRION		50 M		Metal instruments	●
121° C FAST		38 M	3,0	Rubber solids and delicate solids	●
134° C FAST		28 M		Rubber solids and metal solids	●

REFERRING TO KRONOS N18

KRONOS B18	KRONOS S18	KRONOS N18	
fractionated	pre and post vacuum	no vacuum	TYPE OF VACUUM
<p style="text-align: center;">✓</p>			INTEGRATED STEAM GENERATOR
			BUILT-IN PRINTER
			M.T.S MEMORY TEST SYSTEM
			WATER QUALITY SENSOR
			CONNECTION TO EXTERNAL DATA LOGGER
optional			INTERNAL USB DATA LOGGER
505Wx615Dx400H mm			OVERALL MEASURES
ø236x381,5 mm			CHAMBER SIZE
nr. 4 trays: 185Wx285Dx18H mm			TRAY SIZES
230 V ± 10%			SUPPLY VOLTAGE
50-60 Hz			MAINS FREQUENCY
2400 W			RATED POWER

TECHNICAL DATA



KRONOS B18