

# KRONOS 23 L

## RECOMMENDED PROGRAMS

● OPEN ● WRAPPED

B-CLASS

TEMPERATURE AND NAME OF CYCLE	TYPE	TOTAL TIME	WEIGHT	STERILIZABLE MATERIALS AND INSTRUMENTS	WRPD
HELIX TEST - B&D TEST	Test	29 M	0	Steam penetration test	-
VACUUM TEST		15 M		Leak detection test	-
121° C HOLLOW WRAPPED	B	51 M	3,5	Delicate hollow, stainless steel hollow and turbines	●
134° C HOLLOW WRAPPED		49 M		Stainless steel hollow and turbines	●
134° C PRION		60 M		Stainless steel hollow and metal instruments	●
134° C POROUS		58 M	1,0	Porous loads	●
134° C HOLLOW OPEN		41 M	4,5	Stainless steel hollow and turbines	●

REFERRING TO KRONOS B23

S-CLASS

TEMPERATURE AND NAME OF CYCLE	TYPE	TOTAL TIME	WEIGHT	STERILIZABLE MATERIALS AND INSTRUMENTS	WRPD
BOWIE & DICK TEST	Test	35 M	0	Steam penetration test	-
VACUUM TEST		15 M		Leak detection test	-
121° C HOLLOW WRAPPED	S	57 M	2,5	Delicate hollow, stainless steel hollow	●
134° C HOLLOW WRAPPED		55 M		Stainless steel hollow	●
134° C PRION		66 M		Stainless steel hollow and metal instruments	●
134° C POROUS		63 M	0,5	Little porous loads	●
134° C FAST		40 M	3,5	Rubber solids and metal solids	●

REFERRING TO KRONOS S23

N-CLASS

TEMPERATURE AND NAME OF CYCLE	TYPE	TOTAL TIME	WEIGHT	STERILIZABLE MATERIALS AND INSTRUMENTS	WRPD
121° C SOLIDS OPEN	N	58 M	2,5	Rubber solids and delicate solids	●
134° C SOLIDS OPEN		46 M		Rubber solids and metal solids	●
134° C PRION		60 M		Metallic Instruments	●
121° C FAST		48 M	4,5	Rubber solids and metal solids	●
134° C FAST		38 M	Rubber solids and metal solids	●	

REFERRING TO KRONOS N23

KRONOS B23 KRONOS S23 KRONOS N23

fractionated pre and post vacuum no vacuum

TYPE OF VACUUM

INTEGRATED STEAM GENERATOR

BUILT-IN PRINTER

M.T.S MEMORY TEST SYSTEM

WATER QUALITY SENSOR

CONNECTION TO EXTERNAL DATA LOGGER

INTERNAL USB DATA LOGGER

INTEGRATED WATER PURIFICATION SYSTEM "ROSI"

OVERALL MEASURES

CHAMBER SIZE

TRAY SIZES

SUPPLY VOLTAGE

MAINS FREQUENCY

RATED POWER

TECHNICAL DATA



optional

505Wx690Dx400H mm

ø236x530 mm

nr. 4 trays:  
185Wx440Dx19H mm

230 V ± 10%

50-60 Hz

2400 W



KRONOS B23